


multiSUB choice

- Three tray options
- Run up to 210 samples
- Low buffer volumes
- Multichannel pipette compatible combs for speed loading

Typical Applications

Ideal for restriction fragment analysis, sample prep or checking of high numbers of samples



multiSUB Choice Trio includes all 3 tray sizes for optimum versatility and value

The multiSUB Choice offers a wide degree of versatility. Three tray options are available - 15 x 7cm, 15 x 10cm and 15 x 15cm - allowing the choice of one, two or all three gel length options at the time of purchase. Further purchases of additional accessories are no longer required. Maximising comb and tray options allow up to 210 samples to be resolved per gel. The 15cm total run length allows restriction fragment or other close MW sample bands to be easily separated and identified. Speed loading is accomplished using 10, 14, 16, 18, 28 and 30 sample multi-channel pipette compatible combs.

Technical Specifications

Gel dimensions (W x L)	15 x 7cm 15 x 10cm 15 x 15cm
Unit dimensions	26.5 x 17.5 x 9cm (W x L x H)
Max. sample	15 x 7cm tray - 70 samples 15 x 10cm tray - 140 samples 15 x 15cm tray - 210 samples
Buffer volume	500ml
Combs available:	
No. of samples	1, 2, 4, 10, 10MC, 12, 14MC, 16, 18MC, 20, 28MC, 30MC, 35
Thicknesses	0.75, 1, 1.5, 2mm

Ordering Information

MSCHOICE7	multiSUB Choice, 15 x 7cm UV Tray, 2 x 20 sample combs, loading guides and dams		
MSCHOICE10	multiSUB Choice, 15 x 10cm UV Tray, 2 x 20 sample combs, loading guides and dams		
MSCHOICE15	multiSUB Choice, 15 x 15cm UV Tray, 2 x 20 sample combs, loading guides and dams		
MSCHOICETRIO	multiSUB Choice Trio, 15 x 7, 10 & 15cm UV Tray, 2 x 20 sample combs, loading guides and dams		
MS15-UV7	15 x 7cm UV Tray	MS15-NE	Negative Electrode
MS15-UV10	15 x 10cm UV Tray	MS15-LG	Adhesive Loading Guides
MS15-UV15	15 x 15cm UV Tray	MS15-WP	Viewing Platform
MS15-UNVDM	Casting Dams	MS15-UNS	15cm UV Gel Scoop
MS15-PE	Positive Electrode	MSCHOICECP	Cool-Pack and Platform
MS15/20-FC	multiSUB Choice/Maxi Flexi caster		
MSCHOICEBSB	Buffer Saver Blocks p/2, saves 190ml of buffer		
MS15-x-yLG	Combi Comb and Loading Guide, x = comb well no., y = comb thickness		

CODE	DESCRIPTION	SAMPLE VOLUME FOR A 5MM THICK GEL	CODE	DESCRIPTION	SAMPLE VOLUME FOR A 5MM THICK GEL
MS15-1-0.75	Comb Prep 1, Marker 1, 0.75mm thick	371µl	MS15-1-1.5	Comb Prep 1, Marker 1, 1.5mm thick	742µl
MS15-2-0.75	Comb Prep 2, Marker 2, 0.75mm thick	169µl	MS15-2-1.5	Comb Prep 2, Marker 2, 1.5mm thick	338µl
MS15-4-0.75	Comb Prep 4, Marker 2, 0.75mm thick	91µl	MS15-4-1.5	Comb Prep 4, Marker 2, 1.5mm thick	182µl
MS15-10-0.75	Comb 10 sample, 0.75mm thick	24µl	MS15-10-1.5	Comb 10 sample, 1.5mm thick	60µl
MS15-12-0.75	Comb 12 sample, 0.75mm thick	30µl	MS15-12-1.5	Comb 12 sample, 1.5mm thick	61µl
MS15-20-0.75	Comb 20 sample, 0.75mm thick	16µl	MS15-20-1.5	Comb 20 sample, 1.5mm thick	32µl
MS15-35-0.75	Comb 35 sample, 0.75mm thick	7µl	MS15-35-1.5	Comb 35 sample, 1.5mm thick	15µl
MS15-10MC-0.75	Comb 10 sample MC, 0.75mm thick	22µl	MS15-10MC-1.5	Comb 10 sample MC, 1.5mm thick	44µl
MS15-14MC-0.75	Comb 14 sample MC, 0.75mm thick	22µl	MS15-14MC-1.5	Comb 14 sample MC, 1.5mm thick	44µl
MS15-16MC-0.75	Comb 16 sample MC, 0.75mm thick	20µl	MS15-16MC-1.5	Comb 16 sample MC, 1.5mm thick	41µl
MS15-18MC-0.75	Comb 18 sample MC, 0.75mm thick	8µl	MS15-18MC-1.5	Comb 18 sample MC, 1.5mm thick	17µl
MS15-28MC-0.75	Comb 28 sample MC, 0.75mm thick	8µl	MS15-28MC-1.5	Comb 28 sample MC, 1.5mm thick	17µl
MS15-30MC-0.75	Comb 30 sample MC, 0.75mm thick	9µl	MS15-30MC-1.5	Comb 30 sample MC, 1.5mm thick	19µl
MS15-1-1	Comb Prep 1, Marker 1, 1mm thick	495µl	MS15-1-2	Comb Prep 1, Marker 1, 2mm thick	990µl
MS15-2-1	Comb Prep 2, Marker 2, 1mm thick	225µl	MS15-2-2	Comb Prep 2, Marker 2, 2mm thick	450µl
MS15-4-1	Comb Prep 4, Marker 2, 1mm thick	122µl	MS15-4-2	Comb Prep 4, Marker 2, 2mm thick	243µl
MS15-10-1	Comb 10 sample, 1mm thick	45µl	MS15-10-2	Comb 10 sample, 2mm thick	90µl
MS15-12-1	Comb 12 sample, 1mm thick	41µl	MS15-12-2	Comb 12 sample, 2mm thick	81µl
MS15-20-1	Comb 20 sample, 1mm thick	21µl	MS15-20-2	Comb 20 sample, 2mm thick	42µl
MS15-35-1	Comb 35 sample, 1mm thick	10µl	MS15-35-2	Comb 35 sample, 2mm thick	20µl
MS15-10MC-1	Comb 10 sample MC, 1mm thick	29µl	MS15-10MC-2	Comb 10 sample MC, 2mm thick	59µl
MS15-14MC-1	Comb 14 sample MC, 1mm thick	29µl	MS15-14MC-2	Comb 14 sample MC, 2mm thick	59µl
MS15-16MC-1	Comb 16 sample MC, 1mm thick	27µl	MS15-16MC-2	Comb 16 sample MC, 2mm thick	54µl
MS15-18MC-1	Comb 18 sample MC, 1mm thick	11µl	MS15-18MC-2	Comb 18 sample MC, 2mm thick	23µl
MS15-28MC-1	Comb 28 sample MC, 1mm thick	11µl	MS15-28MC-2	Comb 28 sample MC, 2mm thick	23µl
MS15-30MC-1	Comb 30 sample MC, 1mm thick	13µl	MS15-30MC-2	Comb 30 sample MC, 2mm thick	25µl

MC = multichannel pipette compatible